

## Installation Guide for the «Sanatrix» App on Mobile Device

To install the «Sanatrix» app on your mobile device, follow these steps:

1. Ensure that your mobile device has Android version 7 or higher installed.
2. Download the installation file of the app from the same page of the website where this guide is available.
3. Run the downloaded file from the notification window or from the file manager of your mobile device. If your device shows a message that installation from unknown sources is blocked or something similar, proceed to step 4. If the installation occurs, proceed to step 5.
4. Go to the settings of your mobile device and allow installation of unknown apps. Depending on the version of the operating system, this setting may be located along one of the paths listed below or in a similar path:

- Settings > Apps & notifications > Special access > Install unknown apps;
- Settings > Security & privacy > Other permissions > Install unknown apps.

Alternatively, go to app permissions or follow the path below, select the «Package Installer» app, and allow installation from external sources:

- Privacy > Management > Special access > External sources.
5. Next, follow the instructions of the app installation wizard.
  6. Once the app installation is complete, launch it and agree to all permission requests displayed by the program, if you agree. Without agreeing to the requests, the functionality of the app may be limited.

Wishing you successful use of the «Sanatrix» app!